

THESE DIVERSE ATTACHMENTS COME IN THREE STYLES:



Frost – The firm versatile plastic design allows the practitioner easy access of all parts of the anatomy. The Frost is not just for static trigger point application, but has the ability to slide over skin and clothing. It is ideal for long muscles like the IT band, SI issues and forearm flexors.



Cabernet – This firm rubber design works very good for “athletes” and patients that benefit from a firmer technique. The Cabernet is excellent for breaking up local, deeper adhesions and elongating fascia. It’s cylinder design is ideal for the scapular region as well. The attachment absorbs the uncomfortable sting of a firmer therapy technique.



Mint – This soft rubber design works on localized trigger point release. The Mint attachment is best used for static or stationary applications. The soft head is necessary for the “typical” patient and their tolerance to muscle treatment. It provides a big difference in absorbing the discomfort of a firmer head and has the right amount of give to apply the correct pressure for best results.

Percussor Trigger-Point Series

By Erchonia®



Erchonia Corporation
650 Atlantis Rd.
Melbourne, FL 32904



ERCHONIA
World Leaders in Low Level Laser Technology

AS AN INNOVATOR IN CHIROPRACTIC TECHNOLOGY, ERCHONIA IS TAKING OUR PERCUSSOR TECHNOLOGY TO THE NEXT LEVEL

Introducing the New Percussor Trigger-Point Series attachments from Erchonia. They easily attach to any Erchonia Percussor, allowing even the softest tissue technique to become easier and more dynamic. With the Percussor Trigger-Point Series attachments, not only will the practitioner gain the desired result with less time and effort, they make the overall procedure less traumatic for you and your patients. All attachments are designed to be used on a variety of patients of different age, health and fitness status.

For US and International patents and patent pending information go to www.erchonia.com.



Scan with RedLaser



Get the free mobile app for your phone
<http://gettag.mobi>

THE PERCUSSOR TRIGGER-POINT SERIES CAN BE USED TO HELP THE FOLLOWING:

- ▶ Trigger points
- ▶ Trapezius and glute fibers
- ▶ Muscle stripping and break-up adhesions
- ▶ Elongate and release fascia



THE TRIGGER-POINT ADVANTAGE

These three attachments are sold as a set for an affordable price. The Percussor Trigger-Point Series attachments are long lasting, offer incredible value and a strong ROI, helping you put more profit into your practice.

Easy to Use
Storage Stand
(Included with Set)



Be among the first to experience
this innovative new product.

For more information, please
contact Erchonia at

888.242.0571

www.erchonia.com

